

Medorra Adaptive Gown

Adjustable patient gown

A wraparound design with snap closures, an adjustable neckline, and a pass-through telemetry pocket.



Fit and coverage


Waist and neck adjustments accommodate sizes Medium to 3XL. Check fit for the individual wearer.

No waist ties

The back belt wraps to the front and fastens at the selected snap point.

A closer look at the fit

The front-opening gown adjusts at two areas so staff can check coverage and comfort for the individual wearer.



Wraparound waist belt

The belt extends from the back and attaches to snap points at the front waistline. The closure replaces a traditional waist tie.

Adjustable neckline

Multiple neck snap positions allow the opening to be adjusted. Check the selected neck and waist positions for appropriate fit.


Size range: Medium to 3XL. Individual fit depends on body dimensions and the selected configuration.

Features for daily care

Neck adjustment and IV snap sleeves

Snap closures at the neckline and sleeves provide access points in the gown design.


ADJUSTABLE BACK NECK & SHOULDER OPEN FUNCTION:



Pass-through telemetry pocket

The pocket includes an internal pass-through. Metal-free options are available. Confirm the exact configuration with Medorra.

Select the gown according to your facility's care requirements and product instructions. No equipment compatibility is implied.



Plan your gown selection

Talk with Medorra about the configuration, quantities, and ordering options that suit your facility.



Product at a glance

Product	Medorra Adaptive Gown adjustable patient gown
Size range	Medium to 3XL
Closures	Snap waist belt and adjustable neck snaps
Design features	IV snap sleeves and pass-through telemetry pocket
Options	Metal-free configurations. Confirm availability.

Medorra gown inquiries

leon@medorra.com

(845) 263-6158

Explore patient gowns

medorra.com

Medorra Adaptive Gown informational brochure. Confirm current specifications, available configurations and care instructions with Medorra before ordering.